**Silicone-based soft denture liner**  
**Weiche Silikonunterfütterung Für Prothesen**  
**Materiu de revetment souple en silicone pour prothèses dentaires**

# SOFRELINER

REF SR02

**Manufacturer**  
**Tokuyama Dental Corp.**  
 38-9, Taiyū 1-chome, Taiyoku-ku, Tokyo, Japan  
 TEL:+81-3-3835-7201  
 Email:export@tokuyama-dental.co.jp

**EC Rep.**  
**Tokuyama Europe GmbH**  
 Ostrassse 10, 40211 Düsseldorf, Germany  
 Tel: +49-211-175480 Fax: +49-211-357379

**North America**  
**Tokuyama America, Inc.**  
 California  
 TEL:+1(571) 8872  
 TEL:(877)338 3548 (Dental Toll-Free)

## ENGLISH

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## ENGLISH

**SOFRELINER**  
**Read all information, precautions, notes and materials before using the product.**

### PRODUCT DESCRIPTION AND GENERAL INFORMATION

SOFRELINER is an addition-cured silicone chaiside resilient lining material for dentures. Sofreliner has no odor or taste. It is provided as a new automatic cartridge delivery system together with an adhesive primer and other accessories. The material cures in 5 minutes intraorally or in 20 minutes at room temperature (25°C / 77°F).

### CONTENTS

- |   |    |
|---|----|
| 1. Cartridge of paste containing Base (27g) and Catalyst (27g), to be used with the TOKUYAMA DISPENSER II(Sold separately). ..... | x1 |
| 2. Primer (10 mL) .....   | x1 |
| 3. Brushes for Primer .....   | x2 |
| 4. Plastic cap for Primer .....   | x2 |
| 5. Spatula for spreading SOFRELINER after it has been applied to the denture.....   | x1 |
| 6. TD II Mixing tip (S) or (XS) for use with paste cartridge and the TOKUYAMA DISPENSER II .....                                  | x1 |
| x10   |    |
| 7. Steel Bur T 101 .....  | x1 |
| 8. Finishing Wheel .....  | x1 |
| * Note: Only the TOKUYAMA DISPENSER II can be used with SOFRELINER.   |    |

### CLINICAL APPLICATION OF SOFRELINER LABORATORY LINING METHOD

- EXAMINATION OF THE DENTURE
  - If the denture is worn or old, line the denture with TOKUYAMA REBASE II. The TOKUYAMA REBASE II material should be applied in the appropriate thickness so that when the new material is ground down, a newly lined layer of TOKUYAMA REBASE II is left on the denture bearing area allowing SOFRELINER to adhere to it. In this case, make sure that the tissue side of the denture is ground deeply enough, so that the relined SOFRELINER will not be too high on the patients' occlusion after the lining with TOKUYAMA REBASE II.
- FABRICATION OF A GYPSUM MODEL
  - In order to enhance bonding on old or worn dentures, line the denture with TOKUYAMA REBASE II prior to lining with SOFRELINER. In this case, make sure that the tissue side of the denture is ground deeply enough, so that the relining of SOFRELINER will not be too high on the patients' occlusion after the lining with TOKUYAMA REBASE II.

### CLINICAL APPLICATIONS FOR SOFRELINER

- For pain caused by sharp ridges of the alveolar bone and atrophy of the mucosa
- Post retention caused by extreme ridge atrophy
- Major undercuts or maxillary/mandibular tori
- Retention of overdentures before the final attachment is placed

### STORAGE AND DISPOSAL

Keep SOFRELINER in a cool, dark place (0-25°C / 32-77°F) before using. After using SOFRELINER, the mixing tip must be replaced with the cartridge cap, which is attached to the cartridge at the time of purchase. DO NOT use the material if the expiration date has passed. To dispose the paste, dispense unused paste out from the cartridge and harden it before disposal. To dispose the Primer, remove the liquid from the bottle and soak the liquid in tissue paper before disposal.

### CLINICAL PRECAUTIONS

- SOFRELINER is intended for use only by licensed dental care professionals.
- Use SOFRELINER only as instructed herein. DO NOT use SOFRELINER for any other application.
- DO NOT use this product if the safety seals are broken or if the contents have been tampered with.
- Primer contains methylene chloride, which is highly volatile. Avoid long-term exposure, remove the room during and after use, and close the cap tightly after use.
- SOFRELINER is a silicone-based product. DO NOT use SOFRELINER on patients with a history of sensitivity to silicone. Always determine whether your patient is allergic to, or hypersensitive to silicone-based products. If a patient is uncertain whether or not he/she is sensitive to silicone based products, such as SOFRELINER, always perform a sensitivity test on the patient before using this product. Allow a small amount of SOFRELINER paste to set on a dispensing dish. Use this set material as a test patch by placing the material with a sticking plaster on the patient's arm for at least 24 hours, monitoring the test patch area for sensitivity. Advise all patients of the possibility of sensitivity to silicone, and discuss the signs of silicone sensitivity with them. If a patient shows any signs of sensitivity, patient should be instructed to immediately contact their dentist for removal of SOFRELINER from the denture.
- The use of examination gloves (plastic, vinyl) is recommended at all times to avoid the possibility of allergic reactions from SOFRELINER. The use of latex gloves may inhibit the cure of SOFRELINER. If SOFRELINER Primer comes in contact with the examination gloves, immediately dispose of the examination gloves and insert new gloves thoroughly with water.
- AVOID CONTACT WITH EYES. In case of contact with eyes, thoroughly flush eyes with water, and immediately contact a physician.
- Avoid contact of uncured chemical with skin. In case of contact with skin, clean effected area with rubbing alcohol, then thoroughly flush with water.
- Avoid contact of unpolymerized materials with clothing. In case of contact with clothing, clean the material. Allow it to set and peel it off.
- DO NOT allow excess material to flow down the patient's trachea or esophagus. Aspirated or ingested SOFRELINER paste can cause asphyxiation or other injury, even though SOFRELINER paste is non-toxic. Immediately remove any excess lining material from the palatal area to avoid flow of material down the patient's throat. SOFRELINER is not for use on the tongue or in the oral cavity.
- DO NOT mix SOFRELINER with any materials not noted in these instructions.
- Clean the tip of the cartridge with a clean cloth or tissue each time before using a new mixing tip or replacing the cartridge cap.
- The material cures in 5 minutes intraorally or in 20 minutes at room temperature (25°C / 77°F). However, SOFRELINER does not need to be cured until one day after the lining. Therefore, the softness of the material may change slightly during this period. The finishing procedure can be performed immediately after the extraction of the newly lined denture from the mouth or a gypsum model. The newly lined denture can be worn after the finishing procedure.
- Always clean the spatula, brushes and cups with alcohol after each use.
- Instruct the patients in accordance with the PATIENT PRECAUTIONS contained in the last section of these instructions.

### CLINICAL APPLICATION OF SOFRELINER INTRAORAL LINING METHOD

- Remove the denture from the patient's mouth and clean the area to be lined with soap and water. Then air-dry the surface of the denture thoroughly.
- Grind the tissue side of the denture to the required thickness with a bur. Bevel the borders of the denture to allow at least 2 millimeters of material to roll on the borders. Avoid setting the marginal line on the denture border.
- After blowing the ground surface of the denture with strong air stream, roll the Primer into the plastic cup included in the SOFRELINER kit. Using the brush included in the SOFRELINER kit, apply the Primer evenly in a thin coat on the surface of the denture where SOFRELINER is to be lined. Do not fix any spot where the Primer should be applied, since this may cause poor adhesion of SOFRELINER.
- Note: Attach the Drip-cap included in the SOFRELINER kit onto the rim of the Primer bottle to facilitate pouring.
- Load the SOFRELINER paste cartridge into a TOKUYAMA DISPENSER II (sold separately).
- Evenly dispense SOFRELINER paste on the tissue side of the denture. Spread the paste evenly with a spatula across all areas where the Primer was applied.
- Suction any excess saliva from the patient's mouth 1 - 1 1/2 minutes after dispensing SOFRELINER onto the denture. Place the denture in the patient's mouth, and have the patient occlude very gently into centric occlusion to the proper vertical dimension. Check for any excess SOFRELINER in the patient's mouth, and remove immediately if any excess is present. DO NOT allow excess material to flow down the patient's trachea or esophagus. Aspirated or ingested SOFRELINER paste can cause asphyxiation or other injury, even though SOFRELINER paste is non-toxic. If the material flows towards the patient's trachea, immediately wipe excess material from the patient's mouth with a

The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure the product is suitable for an appropriate application before use. Specifications are subjected to change without notification.

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