



Announcement

April 3, 2006

To All Tokuyama Dealers,

Tokuyama America Inc. has introduced the revolutionary flowable composite, **ESTELITE FLOW QUICK**, offering the following benefits:

- **Super quick curing:** 5-10 sec. with any curing light
- **Outstanding esthetics:** submicron spherical filler gives excellent polishability and long lasting gloss retention.
- **Superior mechanical strengths:** designed for both anterior and posterior restorations
- **Low shrinkage:** high filler content (71%/weight) offers one of the lowest shrinkage percentage.
- **Chameleon effect:** material's optimal translucency and submicron spherical filler make it easy to match and blend.
- **Minimum color change:** controlled refractive index of filler minimizes the color change of material after curing.

ESTELITE FLOW QUICK is also a radiopaque and available in 17 shades.

Please contact **Tokuyama America Inc.** for any inquiries concerning Tokuyama Dental products.

Tokuyama America Inc.
1799 Old Bayshore Hwy, Suite 168,
Burlingame, CA 94010-1319
Tel: 1-877-378-3548 or (650) 571-8872
Fax: (650) 571-8037
info@tokuyama-us.com