



LUXURY RESTORATIONS

Introducing the next generation *ESTELITE Sigma Quick...* PALFIQUE LX5. Delivering lifelike aesthetic restorations, infused with our Supra-Nano spherical fillers, that was engineered by Tokuyama Dental. Tokuyama Dental's composite resin is recognized world-wide, and features outstanding polishability, wide shade matching range and high wear resistance. PALFIQUE LX5 can meet everything in your practice.

RESIN-BASED DENTAL RESTORATIVE MATERIAL

PALFIQUE LX5




BENEFITS

- Quick curing time - 10 sec. with a halogen light ($\geq 400\text{mW/cm}^2$)
- Extended working time - 90 sec. under ambient light (10,000 lx)
- Outstanding polishability
- High gloss retention over time
- Wide shade matching range (chameleon effect)
- High wear resistance
- Less wear to opposing teeth
- Low shrinkage
- Good radiopacity
- Available in 21 shades

CLINICAL CASES

CLASS II



Before




After

Case courtesy of Hirofumi Tashiro, D.D.S., Ph.D.

CLASS IV




Before



After

Case courtesy of Tomohiro Takagaki, D.D.S., Ph.D.

EXCELLENT PHYSICAL PROPERTIES




** Source: Tokuyama Dental R&D
* Not a registered trademark of Tokuyama Dental


SUPERIOR & LONG LASTING GLOSS

The Supra-Nano Spherical fillers in PALFIQUE LX5 provide a very smooth surface for a superior shine, and high gloss retention.

FILLER SEM IMAGES




PALFIQUE LX5




Venus Pearl

GLOSS RETENTION



PALFIQUE LX5



Venus Pearl

PRODUCT PACKAGES

Please check with your local distributor for available shades and packagings.

PALFIQUE LX5 **Kit**

PALFIQUE LX5 SYRINGE SYSTEM KIT

Contents: 6 PALFIQUE LX5 2.0mL (3.8g) syringes (A1, A2, A3, A3.5, OA3, OPA2),
3 ESTELITE FLOW QUICK 1.0mL (1.8g) syringes (A2, A3, OA3), 27 ESTELITE FLOW QUICK Syringe Tips,
1 PALFIQUE BOND (5mL), 1 Dispensing Well, 50 Disposable Applicator (Superfine), 1 shade guide

PALFIQUE LX5 SYRINGE INTRODUCTORY KIT

Contents: 6 PALFIQUE LX5 2.0mL (3.8g) syringes (A1, A2, A3, A3.5, OA2, OPA2),
1 PALFIQUE BOND (1mL), 1 Dispensing Well, 25 Disposable Applicator (Superfine), 1 shade guide

PALFIQUE LX5 3 SYRINGE KIT

Contents: 3 PALFIQUE LX5 2.0mL (3.8g) syringes (A2, A3, OA3)

PALFIQUE LX5 **Refill**

PALFIQUE LX5 SYRINGE REFILL

Contents: 1 syringe 2.0mL (3.8g)
Available shades: A1, A2, A3, A3.5, A4, A5, B1, B2, B3, B4, C1, C2, C3, OA1, OA2, OA3, OPA2, BW (Bleach White), WE (White Enamel),
CE (Clear Enamel), SW (Snow White)

PALFIQUE LX5 PLT REFILL

Contents: 20 PLTs (0.2g each)
Available shades: A1, A2, A3, A3.5, A4, A5, B1, B2, B3, B4, C1, C2, C3, OA1, OA2, OA3, OPA2, BW (Bleach White), WE (White Enamel),
CE (Clear Enamel), SW (Snow White)

Specifications are subject to change without notice based on Manufacturers discretion.

Manufactured by
Tokuyama Dental Corporation

38-9, Taitou 1-chome, Taitou-ku, Tokyo 110-0016, Japan
TEL: +81-3-3835-2261 FAX: +81-3-3835-2265
Email: tokuyama-inquiry@tokuyama-dental.co.jp
URL: www.tokuyama-dental.com